

Biological Microscope

Model: STM-2030FB/FT

Instruction Manual



This instruction manual is for STM-2030FB/FT biological microscope. To insure safety and obtain optimum performance and familiarize yourself fully with the use of this microscope. We recommend that you read the manual thoroughly before operating the microscope, Attain this manual instruction in an easily accessible place near the microscope for the further reference.

Microscope Optical Part

CONTENTS

User Notice.....	2
1.Components	4
2.Assembly	5
2.1 Assembly Diagram.....	5
2.2 Assembly Steps.....	6
3. Adjustment & Operation	10
3.1 Adjustment Sets	10
3.2 Operation	11
4. Specification Table.....	14
5. Outfit	15
6. Troubleshooting Guide	16

User Notice

I . Safety Notes


1. Carefully open the box, avoid the accessories, like lens, dropping to ground and being damaged.
2. Do keep the instrument out of direct sunlight, high temperature or humidity, dusty and easy shaking environment. Make sure the stage is smooth, horizontal and firm enough.
3. When moving the instrument, please use two hands to grip with the two sides of the microscope body.
4. When running, the lamp house and nearby parts will be very hot. Please ensure there is enough cooling room for them.
5. Make sure the instrument is earthed, to avoid lighting strike.
6. For safety, be sure the main switch is in "O"(off) state before replacing the halogen/LED lamp or the fuse, then cut off the power, and do the operation after the lamp bulb and the lamp house completely cool down. **(Specified: Halogen Lamp6V/20W)**
7. Check the input voltage: be sure the input voltage which signed in the back of the microscope is consistent with the power supply voltage, or it will bring a serious damage to the instrument.
8. Use the factory supplied power cord, please.

II. Maintenance

1. All the lenses have been well checked and adjusted. It is forbidden to disassemble them yourself.
2. The nosepiece and coarse/fine focus unit have a compact and precise frame; please don't disassemble them as possibly as you can.
3. The nosepiece and coarse/fine focus unit have a compact and precise frame; please don't disassemble them as possibly as you can.
4. Keep the instrument clean, wipe dust regularly, and be attention to avoid contaminating the optical elements especially.
5. The contaminations on the prism, as finger mark and oil, could be gently wiped with a piece of soft cloth or tissue paper, gauze which has been immersed in pure alcohol or aether. **(Note that the alcohol and ether are highly flammable, do keep them away from the fire or potential sources of electrical sparks, and use them in a drafty room as possible as you can.)**
6. Do not attempt to use organic solvents to clean the microscope components other than the glass components. To clean them, use a lint-free, soft cloth slightly moistened with a diluted neutral detergent.

7. When using, if the microscope is splashed by liquid, cut off the power at once, and wipe up the moisture.
8. Do not disassemble any parts of the microscope, which will affect the function or decline the performance of the microscope.
9. Place the instrument in a cool, dry position. When not using the microscope, keep it covered with a dust cover. Make sure the lamp socket is cool before covering the microscope.

1.Components




STM-2030B Biological Microscope

2.Assembly

2.1 Assembly Diagram

The following figure shows the installation sequence of the components. The number in the figure shows the assembly steps.

- ★ Before installing, be sure every components is clean, do not score any parts or glass surface.
- ★ Keep well with hexagon wrench provided. When replacing the components, you will need it again.



2.2 Assembly Steps




Fig.1




Fig.2




Fig.3




Fig.4

2.2.1 Installing Binocular viewing head (Fig.1, 2)

Insert the digital viewing head into the microscope head, turn into the right position, then screw down the bolt ① to fix it.

2.2.2 Installing the eyepieces (Fig.3, Fig.4)

Insert the eyepieces into the eyepiece tube until they are against each other as shown in Fig.4.

NOTE:

Operation Conditions:

1. Temperature: 0°C ~ 40°C, Maximum Relative Humidity: 85%.
2. High Temperature: High Temperature and humidity will result in a mildewing, dew and even ruinous instrument.
3. Avoid placing the instrument in a dusty environment. When ending your microscope operation, please cover it with the dust cap.
4. Please lay the microscope in a plan and stable position.




Fig.5




Fig.6




Fig.7

2.2.3 Installing objectives (Fig.5& 6)

1. Adjusting the coarse focus knob until the support device of the mechanical stage reaches its low limit position.
2. Screw the lowest magnification objective into the nosepiece from the left or the right side, then revolve the nosepiece clockwise and mount other objectives by the sequence of low to high magnification

✧ Installing objective this way will make the change of magnification to be easier during using.

- ★ Clean the objective regularly, for lens is susceptible to dust.
- ★ When operating, use 10×magnification objective to search and focus specimen firstly, then replace with higher magnification objective if necessary.
- ★ When replacing the objective, slowly turn the nosepiece until you hear “clicked”, which means the objective is in the required position--the light path center.

2.2.4 Installing the color filters (Fig.7)

1. Turn the condenser bracket① out at the direction of arrow in Fig.7
2. Put the required filters② into the holder on the bracket, and then turn the bracket back to the right position.

★ Baby blue and green filters are available in standard outfit.




Fig.8




Fig.9




Fig.10

2.2.5 Connecting the power cord (Fig.8, 9, 10)

★ The cable and cords are vulnerable when bent or twisted, never subject the power cord to excessive force.

1. Turn the main switch① to “O” (off) state before connecting the power cord.
2. Insert the power plugs ② into the power jack③ of the microscope; make sure the connection is well.
3. Plug the power cord ④ into the power supply receptacle safely. Make sure the connection is well.

★ Do use the supplied power cord all the time. If lost or damaged, please select the same standard cord.

2.2.6 Replacing the Fuse (Fig.8, 9, 10)

Do remember to set the main switch① to the state of “O” (OFF) and unplug the power cord ② before replacing the fuse. Rotate the fuse kits⑤ out of the holder, replace with a new fuse, then rotate it back to the holder again.

★ The rating of the new fuse should be the same as the old one.




Fig.11




Fig.12

2.2.7 Installing and replacing the lamp (Fig.11, 12)

✧ Please use the specified halogen Lamp (6V20W) or LED lamp. Replace the lamp as follows if necessary.

1. Please set the main switch to “O” (off) state before replacing, and make sure the bulb, the lamp room and periphery are all cool enough to carry no burn. Then, you can do your replacing.
2. Loose the bolt ① and open the window② on the bottom of the microscope base with “—”type screwdriver.
3. Pull out the old bulb③ , hold the new bulb after you wrap it with gauze or other protection materials and insert its pin as deeply as possible into the jack in the lamp holder.
4. Close the window and tighten the bolt ① .

★ Please insert the bulb gently, or it will be damaged by excessive extrusion.

★ Do not touch the halogen bulb with bare hands. It will shorten the service life or cause it to burst. If you leave fingerprints on the surface carelessly, clean it with a piece of dry soft cloth.

3. Adjustment & Operation

3.1 Adjustment Sets (Fig.14, Fig.15)




Fig.14




Fig.15

3.2 Operation




Fig.16




Fig.17




Fig.18



Fig.19

3.2.1 Adjusting the brightness (Fig.16)

1. Connect the power, turn on the main switch ① (shown in the figure) which on the bottom side of the base to “—”(on).
2. Turning the brightness adjustment knob ② clockwise, the voltage decline, and the brightness weaken; Whereas turning at the opposite direction, the voltage raise, and the brightness strengthen.

★ Using the microscope at a lower voltage can prolong the service life of the bulb.

3.2.2 Placing the specimen (Fig.17)

1. Place the specimen ③ on the center of the stage, and then nip it with the specimen holder ④.
2. Turn the transversal and longitudinal adjustment knobs which on the mechanical ruler to move the specimen onto the required position.

★ Be careful when changing the objective. If you finish the observation with the short working distance objective, and want to change another one, be careful of not letting the objective touch the specimen.

3.3.3 Focusing the specimen (Fig.18, 19)

1. Focus the specimen with 10X objective. To avoid the objective touching the specimen during focusing, you should raise the mechanical stage to let the specimen close to the objective at first, then slowly separate them to bring the specimen to focus.
2. Turn the coarse focus knob ① conversely to lower the specimen and search images in the 10×ocular simultaneously, and then use the fine knob ② to make focus. After that, you can replace with other magnification objectives safely, and focus without the risk of damaging the specimen.
3. When the specimen is focused, tighten the locking set (in Fig.14) to give a height limit to the stage. In this way, the stage only can lower, avoiding false operation and damage to the specimen.




Fig.20

- The tight tension of the coarse focus knob has already been adjusted before leaving factory. If loosen (e.g. the stage slip down by its weight), please screw the intention adjustment collar③ to the right position by the supplied spanner.

3.3.4 Condenser Adjustment (Fig.20)

Turn the condenser focus knob to move the condenser up and down. Raise the condenser when using the high magnification objective, and descend it when using the low magnification one.

- ★ The condenser and the objective are coaxial. It has been adjusted before leaving factory, so the user needn't to adjust them by self (the distance between the top of the condenser and the stage should be in the range of 0.03mm~0.4mm.)
- ★ The highest position of the condenser has been adjusted too. It also needn't any user's operation.




Fig.21

3.3.5 Aperture Iris Diaphragm Adjustment (Fig.20, 21)

Turn the aperture iris diaphragm lever ① to adjust the aperture iris diaphragm.

Generally, setting the aperture iris diaphragm to 70-80% of the N.A. of the objective in use will provide an image with good contrast.

- If the size of the aperture diaphragm minified, the brightness and the resolution declined, while the contrast and the depth of field increased; In other words, if the size largen, the brightness and the resolution improved, but the contrast and the depth of field declined.
- Generally, setting the size of the condenser aperture diaphragm at 70%~80% of the numerical aperture, you can obtain a clear image with enough contrast. If the open of the aperture diaphragm is too small, the resolution were very low, so please don't adjust it too small.

- Minify the aperture below 60% of the objective's numerical aperture unless in a special case, for instance, observing an almost transparent specimen.
- The numerical aperture is marked on the objective. For example, the mark " 10/0.25 " means the magnification is 10×, and the numerical aperture is 0.25.
- If you want to observe the image of the aperture iris diaphragm, remove one eyepiece and look through the tube. You will see a dark circle encroaching on the bottom of the tube.




Fig.22

3.3.6 Adjusting the Interpupillary Distance (Fig.22)

The interpupillary distance range:

55mm~75mm.

While looking through the eyepieces, move both sliders to adjust interpupillary distance until the left and right fields of view coincide completely. The value on the right slider is your interpupillary distance.

3.3.7 Adjusting the Diopter (Fig.23)

Turn the diopter rings② on both eyepieces to align the number on the rings (which is corresponding to your interpupillary distance) with the white dot ③ on both sliders to bring the specimen into focus.




Fig.23

★ The diopter range of the eyepiece is ± 5 diopter. The number aligned to the line on the viewing head is the diopter in use.

4. Specification Table

4.1 Main specifications

Mechanical Tube Length	160mm
Viewing Head	Sliding Monocular/Binocular/Trinocular Head, 45°Inclined, Interpupillary Distance 55-75 mm
Eyepiece	Field of view: ϕ 18mm
Nosepiece	External Quadruple Nosepiece
Objective	Achromatic objectives 4 \times , 10 \times , 40 \times , 100 \times (immersion oil)
Focusing	Coaxial Coarse & Fine Adjustment, Moving Range 30mm, Fine Division 0.002mm.
Condenser	Abbe Condenser, NA=1.25
Stage	Double Layers Mechanical Stage 140mm \times 140mm, Moving Range 75 \times 50mm
Illumination	LED illumination or 6V/20W Halogen Lamp

4.2 Eyepieces and Objectives

1. Achromatic Objectives

Magnification	Numerical Aperture (NA)	Thickness of glass slide (mm)	Focal length (mm)	Working Distance(mm)	Type
4 \times	0.10	0.17	31.05	18	Dry
10 \times	0.25	0.17	17.13	6.5	Dry
40 \times	0.65	0.17	4.65	0.53	Dry
100 \times	1.25	0.17	2.906	0.13	Oil

Plan Achromatic Objectives are optional

Designation	Magnification	Numerical Aperture	Working Distance
Plan Achromatic Objective	4X	0.10	15.8mm
	10X	0.25	12.2 mm
	40X(S)	0.65	0.37 mm
	100X(Oil,S)	1.25	0.13 mm

1. Eyepieces

Designation	Magnification	Field of view(mm)	Focal Length(mm)
Wide Field	10X	Φ18	24.94mm
Plan Field	16X	Φ11	15.58mm

4.3 Total Magnification

Eyepiece	10×	10×	10×	10×
Objective	4×	10×	40×	100×
Total Magnification	40×	100×	400×	1000×

5. Outfit

Component Name	Specification	Quantity	Standard Outfit
Main body	Main Standard	1	Optional
	Double Layers Mechanical Stage	1	Optional
	Condenser Holder	1	Optional
Viewing Head	Digital binocular head	1	Optional
Condenser	Abbe Condenser, NA=1.25	1	Optional
Nosepiece	Quadruple	1	Optional
Illumination	LED illumination	1	Optional
	6V20W Halogen Lamp	1	Optional
	Spare lamp (6V20W Halogen lamp)	2	Optional
	Spare fuse(2A)	2	Optional
Eyepieces	10×Plan Eyepieces	2	Optional
Objectives	Achromatic objective 4×	1	Optional
	Achromatic objective 10×	1	Optional
	Achromatic objective 40×	1	Optional
	Achromatic objective 100×(oil, spring)	1	Optional
Condenser	Bright Field Condenser with Adjustable Iris Diaphragm	1	Optional
Filter	Blue, Green	1 ea.	Optional

6. Troubleshooting Guide

1. Optical system

TROUBLE	CAUSE	SOLUTION
The edge of the field of view is dark or the brightness is not uniform	The nosepiece is not in the located position (objective and light path not coaxial)	Locate the nosepiece properly where it clicks
	The image of filament is not centered	Center the filament
	A lens (the objective, condenser, eyepiece or collector) is dirty.	Clean it thoroughly
Find dust and stain in the field of view	There are stains on the lens (including condenser, objective, eyepiece and collector)	Clean it up
	There are stains on the specimen	Clean it up
	The position of the condenser is too low	Loosen the condenser's locking bolt, adjust the condenser to the right position
The image is defocused (low resolution \ contrast)	There is no cover slip on the specimen	Add coverslip
	The cover slip is too thick or too thin	Use the standard coverslip (0.17mm)
	The specimen is placed inversely	Reversal it back
	There was oil on the dry objective(easily happened in 40X objective)	Clean it up
	There are stains on the lens (including condenser, objective, eyepiece and collector)	Clean it up
	didn't use oil for the oil objective	Use immerse oil
	There was bleb in the oil	Eliminate the bleb
	Use a unsuitable oil	Change to the specified one
	The size of the aperture diaphragm is too big	Minify it
	There are stains on the incident lens of the binocular tube	Clean it up
	The size of the aperture diaphragm is too small	Open it up
The position of the condenser is too low	Adjust the position	
One side of the image is dark	The condenser is not in the center of the field of view\the condenser inclines	Install the condenser again and adjust the center carefully by centering the bolt
	The nosepiece is not in the right position	Turning it until it reach the "clicked" position
	The specimen is floating	Fix it
The image shift during focusing	The specimen slips on the stage	Fix it
	The nosepiece is not in the right position	Turn it to the " clicked "position

The image is a little yellow	Not use the blue color filter	Use the blue filter
The brightness is not enough	The size of the aperture diaphragm is too small	Adjust again
	The position of the condenser is too low	Adjust the position
	There are stains on the lens (including condenser, objective, eyepiece and collector)	Clean it up

2. Mechanical system

TROUBLE	CAUSE	SOLUTION
The image can not focus when using high magnification objective	The specimen is placed inversely The coverslip is too thick	Turn inversely Use the standard coverslip (0.17 mm)
The objective touch the specimen when changed from low magnification to the higher magnification	The specimen is placed inversely The coverslip is too thick	Turn inversely Use the standard coverslip(0.17 mm)
The specimen is not easy to move	The specimen holder is not fixed	Fix it
The binocular image is not coincident	The interpupillary distance is not correct	Adjust it
Eyes are too tired	No diopter adjustment	Adjust the diopter correctly
	The brightness is not suitable	Adjust the voltage of the lamp

3. Electrical system

TROUBLE	CAUSE	SOLUTION
The lamp can't light when the switch is turned on	No power	Check the connection of the power cord
	The bulb is not inserted	Insert it correctly
	The bulb burns out	Replace it
The lamp burns out suddenly	Use a substandard lamp The voltage is too high	Use the specified lamp to replace, if the problem is not solved, contact with the service department
The brightness is not enough	Use a substandard lamp The voltage is too low	Use the specified lamp increase the voltage
The bulb flickers or the brightness is vertiginous	The bulb is going to burn out	Replace it
	The bulb is not entirely inserted into the holder	Check and insert it again

Epi-Fluorescent Attachment Part

Contents

	Contents	2
	User Notices	3
1	Components Name	5
2	Installation	7
	2.1 preparation	7
	2.2 Mounting the Mercury Burner	7
	2.3 Mounting Protection barrier	7
	2.4 Connecting Cords	8
	2.5 Fuse Replacement	8
3	Adjustment & Operation	9
	3.1 Name of Components	9
	3.2 Operation	11
4	Technical Specifications	13

User Notices






The epi-fluorescent attachment is designed for use with our biological Microscope such as STM-2030B/T and STM-2030B/T.

Safety Note

1. The epi-fluorescent attachment is a precision instrument. Carefully open the box, avoid the accessories dropping to ground and being damaged.
2. Do keep the instrument out of direct sunlight, high temperature or humidity, dusty and easy shaking environment.
3. Make certain that the burner is installed firmly and that all cords are connected firmly.
4. Do not open the lamp housing while it is turned on or for at least 10 minutes after it has been turned off. Lamp housing parts are extremely hot and cause burns if touched.
5. Always be sure to ground (earth) the equipment.
6. Verify that the voltage and the frequency of the AC mains outlet match the setting of the voltage switch and the frequency switch on the rear of the power supply units.
7. Always use the power cord provided and make sure that the main switch is moved to "O"(OFF) before connecting the power cord plug to the wall outlet.
8. To prevent any hazard, always turn the main switch on the power supply unit to "O" (OFF), unplug the power cord plug from the mains outlet before replacing the burner or the fuse, and wait for at least 10 minutes before replacing the burner. (be sure to use a GCQ-100 burner.)
9. To prevent obstruction of the air flow, it is important to make sure to leave enough space around and above the lamp housing.


Safety Symbols

The following symbols are found on the system. Study the meaning of the symbols and always use the equipment in the safest possible manner.

Symbol	Explanation
	Indicates that the surface becomes hot, and should not be touched with bare hands.
	Indicates that the high voltage (upper 1KV) inside, improper handling could result in an electric shock to the use.
	Before use, carefully read the user manual. Improper handling could result in personal injury to the user and/or damage to the equipment.
	Indicates that the main switch is ON.
	Indicates that the main switch is OFF.

- This manual is written for the epi-fluorescent attachment and before equipping it with other ordinary microscope, be sure to learn how to use the microscope.

Maintenance and Storage

1. Clean all glass components by wiping gently with gauze. To remove fingerprints or oil smudges, wipe with gauze slightly moistened with a mixture of ether (70%) and alcohol (30%).
 Since solvents such as ether and alcohol are highly flammable, they must be handled carefully. Be sure to keep these chemicals away from open flames or potential sources of electrical sparks—for example, electrical equipment that is being switched on or off. Also remember to always use these chemicals only in a well-ventilated room.
2. Do not attempt to use organic solvents to clean the non-optical component of the equipment. To clean these, use a lint-free, soft cloth lightly moistened with a diluted neutral detergent.
3. Do not disassemble any part of the power supply unit as malfunction or damage may occur.
4. In order not to impair the safety of the equipment, replace the burner when the counter of NFP-1 indicates “100.00” hours. To prevent any hazard, always turn the main switch on the power supply unit to “O” (OFF), unplug the power cord plug from the mains outlet, and wait for at least 10 minutes before replacing the burner. High-pressure gas is sealed within the mercury burner. Thus, if it is continued to be used past its service life expectancy, the glass tube may deform and may sometimes rupture.

1. Components Name

●Epi-fluorescent Attachment: (Fig.1)

- ① Main body of the Epi-fluorescent Attachment
- ② Power supply unit
- ③ Power cord (please use the power cord provided)
- ④ A GCQ-100 mercury burner
- ⑤ Fuses (8A)
- ⑥ Fluorescence free objective 4×, 10×, 40×, 100×
- ⑦ Clamping screw





Fig.1

Assembly Diagram:

● Fluorescent Microscope:



Assembly of Fluorescent Microscope:

1. Loosen the setscrew ① and take the Trinocular Viewing Head ③ from the body of biological microscope.
2. Insert the epi-fluorescent attachment into the biological microscope correctly and tighten the setscrew ① until it is installed firmly.
3. Insert the Trinocular Viewing Head ③ into the epi-fluorescent attachment correctly and tighten the setscrew ② until it is installed firmly.
4. Replace the old objectives with the new fluorescence free objectives




Fig.2




Fig.3




Fig.4

2.1 Preparation

Carefully open the box, remove all packing material and take the attachment out.

2.2 Mounting the Mercury Burner

(Fig.2 and Fig.3)

1. Loosen the burner socket clamping screw ①, and remove the burner socket. (fig.3)
2. After removing the foam backstop②, Securely attach the + pole of the specified mercury burner③ to the lower mount and the - pole to the upper mount, then tighten the socket clamping screws④.
3. Close the burner socket with burner④ into the original position and tighten the socket clamping screw①.

- Be sure to use a GCQ-100 burner.
- Never subject the burner to excessive force when mounting the Mercury Burner.
- Be careful and avoid leaving fingerprints or dirt on the mercury burner. Attached stain may cause distortion in glass which could result in a ruptured burner. If stained, clean by wiping gently with gauze slightly moistened with a mixture of ether (70%) and alcohol (30%).

★ To prevent any hazard, always turn the main switch on the power supply unit to “O” (OFF), unplug the power cord plug from the mains outlet, and wait for at least 10 minutes before replacing the burner.

2.3 Mounting Protection Barrier (Fig.4)

Install the protection barrier⑥on the attachment by tightening the clamping screw⑤.




Fig.5

2.4 Cable and Cord Connections (Fig.5)

1. make sure that the main switch ① of the power supply is set to “O” (OFF) before connecting cables.
2. Plug the connector ② from the burner socket securely into the connector on the power supply unit and make sure the cord is correctly connected.
3. Connect the power cord connector ③ into connector on the power supply unit and make sure the cord is correctly connected.

- Verify that the voltage and the frequency of the AC mains outlet match the setting of the voltage switch and the frequency switch on the rear of the power supply units and improper setting may degrade burner performance , or in the worst case(although very rare), cause the burner to explode.

- It is better to use the power cord provided and the same type power cord should be used if you lose or damage the old one.




Fig.6


2.5 Fuse Replacement (Fig.6)

- 1 Set the main switch ① to “O” (OFF) and unplug the power cord before replacing fuses.
- 2 Using a flat-blade screwdriver, Remove each of the fuses holders ④ by turning it counter- clockwise and pulling out.
- 3 Replace both fuses with new ones.


- Always use the designated fuses (8A). And make sure the voltage of the fuse match the voltage of the AC mains outlet.



3. Adjustment & Operation

3.1 Name of Components



- The reflected light fluorescent mirrors for B-excitation and G-excitation have been installed in the filter system at the factory.
- The assistant excitation filters and cutoff filters for fluorescence are choose and buy.





3.2 Operation3

3.2.1 Preparation

1. Verify that the voltage and the frequency of the AC mains outlet match the setting of the voltage switch and the frequency switch on the rear of the power supply units.
2. Make sure the cord is connected firmly.
3. When transmitted light observation is required, pull out the filter system and make the hole in the light path.
4. Be sure to use immersion oil when using fluorescent free objectives.
5. When it is required to interrupt observation for a short period, use the shield in the assistant excitation filter part. (Repeated on-off of the mercury burner will shorten its service life considerably)
6. Precautions on the specimen color fading:

The system employs high-intensity excitation light to enable bright observation of dark fluorescent specimens. As a result, if high-power objectives are used frequently color fading of the specimen occurs early, degrading the view (contrast) of fluorescent images. So it is effective to use the shutter frequently to avoid illuminating the specimen for a longer period than required.

Color fading of the specimen can also be delayed using commercially available color fading preventing agent (DABCO, etc). The use of color fading preventing agent is recommended when you perform high-magnification observation frequently.

★ Note that color fading preventing agent cannot be used with certain specimens

3.2.2 Select Fluorescent Filter Combination


Select fluorescent filters combination according to the fluorescent dye you use.

excitation	applications
U	·Auto-fluorescence observation ·DAPI: DNA ·Hoechst 332528, 33342: Chromosome
V	·Catecholamines ·5-hydroxy tryptamine ·Tetracycline: Skeleton, Teeth
B	·FITC: Fluorescent antibody method ·Acidine orange: DNA, RNA ·Auramine: Tubercle bacillus ·EGFP, S65T, RSGFP
G	·Rhodamine, TRITC: Fluorescent antibody method ·Propidium iodide: DNA ·RFP

3.2.3 Switch on Electrical Source

Set the main switch of the power supply unit to "1" (ON). The arc will stabilize in 5 to 10 minutes after ignition.

- Some mercury burners may not ignite the first time the power is turned ON due to variance in production, and the safety mechanism in the starter in such a case. If this occurs, set the main switch to “1” (ON), then press the starter reset switch on the front panel of the power supply and between 1 to 4 seconds are required for igniting the burner. Repeat as necessary.
- To avoid shortening the burner life, do not turn the burner off within 15 minutes after ignition.
- The burner cannot be re-ignited for about 10 minutes, that is, until the mercury vapor inside it has cooled down and condensed to liquid.
- Ensure that the hour counter is reset to “000.00” after replacement of the burner. And you can insert a thin object such as a mechanical pencil tip into the reset hole on the front panel of the power supply unit to press the internal switch.



3.2.4 Adjustment the Optics System (Fig.1)

1. Adjust the filter system and put the protection filter in the light path. (To avoiding burning the eye, strictly prohibit from putting other hole in the light path when adjusting the optics system)
2. Put the hole which is in the center of the assistant excitation filter part in the light path.
3. According to the image on the barrier for focusing and centering, make the burner in the center of the light by adjusting the burner centering screw and the collector focusing knob.
4. When it is required to interrupt observation for a short period, use the shield in the assistant excitation filter part.

4. Technical Specifications:

Epi-Fluorescent Illumination	Fluorescent Filter Specifications	
	B Exciting Light Filter System	Standard
	G Exciting Light Filter System	Standard
	U Exciting Light Filter System	Optional
	V Exciting Light Filter System	Optional
Illumination	100WHBO Ultra Hi-voltage Spherical Mercury Lamp	Standard
Protection Barrier	Barrier to Resist the Ultraviolet Light	Standard
Power	Power supplier NFP-1, AC Input 220V/110V (Interchangeable), Digital Display and Timer	Standard
	Power Supplier NFP-A, AC Input 220V, Indicator Display	
Objective	Fluorescence Free Objective 4X	Standard
	Fluorescence Free Objective 10X	Standard
	Fluorescence Free Objective 40X	Standard
	Fluorescence Free Objective 100X	Standard
Immersion Oil	Fluorescence Free Oil	Standard
Vertical Illumination	StandardAchromatic Optics System StandardFiltering System	Standard B and G Exciting Light Filter System
	StandardObservation Methods ① Fluorescence ② Transmitted Light	
Mercury Lamp Housing	StandardMercury lamp housing 100w StandardMercury Burner GCQ100	Standard Standard
Operating Environment	StandardIndoor Use StandardAttitude: Max. 2000m StandardAmbient Temperature: 5°Cto40°C(41°Fto104°F) StandardMaximum Relative Humidity 80% for Temperature Up to 31°C(88°F) Decreasing linearly through 70% at 34°C (93°F),60 %at 37°C(99°F) to 50% relative humidity at 40°C(104°F) StandardMain supply voltage fluctuations not to exceed ±10% of the nominal voltage StandardPollution Degree 2(in accordance with IEC 664) StandardInstallation/ Over voltage Category II(in accordance with IEC 664)	