

# STALWART

## BOD Analyzer SBD-806



## Features

- Use international mercury free pressure difference method to test BOD content.
- 6 independent testing channels, with no interference between cultivation bottles.
- Real time display of testing results during testing.
- Automatically complete

## Specification

Model	SBD-806
Test Range	0-4000 mg/L
Test Error	It meets the accuracy testing standard for water quality BOD5 (glucose glutamic acid standard solution BOD5=180-230mg/L)
Test Duration	1-30 days
Test Quantity	≤ 6 pieces
Recording Interval	24 minutes -12 hours/time
Storage quantity	Detailed data: Save the BOD values and biochemical reaction curves of the latest set of analyticals samples.
	Historical data: Can save 20 sets of final test data.
Culture flask solvent	580ml
Culture temperature	20 ±1°C
Print	Wireless print
Power	110/220V, 50/60Hz